

Anatomico and Comfort
Head support systems

capitis Confort

The ideal comfort for the head

Pads



SMALL Pad
PMDC04P



MEDIUM Pad
PMDC04M



LARGE Pad
PMDC04G



Anti-decubitus
padding



Removable cover
lycra/neoprene air-x-change



Breathable



Fire retardant



Multi-adjustable
and pivoting

Hardware



LONG Horizontal Bar





SHORT Horizontal Bar
PMDC02C



Vertical Bar

Mounting examples



Configurations

Capitis confort including hardware with long horizontal bar and detachable mounting unit



ART.	SIZE
PMDC01P	Small
PMDC01M	Medium
PMDC01G	Large

Capitis confort including hardware with short horizontal bar and detachable mounting unit

ART.	SIZE
	Small
PMDC01MC	Medium
PMDC01GC	Large



Technical features

CAPITIS CONFORT MEASUREMENTS expressed in centimeters



ART	A	B	C	D	E	F
PMDC04P	21	14	13	15	08	04
PMDC04M	30	14	13	24	10	05
PMDC04G	42	16	14	34	14	07

CAPITIS ANATOMICO MEASUREMENTS expressed in centimeters



ART	A	B	C	D	E	F
PMDC03XP	18	13	08	13	04	02
PMDC03P	24	17	12	19	06	04
PMDC03M	28	19	13	21	08	06
PMDC03G	32	25	16	23	08	08

"For partially sighted users this document can be visualized in PDF format on the web site www.promedicare.eu with the option to enlarge the text."

capitis Anatomico

Contour and support to the movement

Pads



EXTRA SMALL Pad



SMALL Pad



MEDIUM Pad



LARGE Pad



Low Profile and
Soft padding



Removable Cover
lycra/neoprene air-x-change



Breathable



Fire retardant



Multi-adjustable
e pivoting




Detachable
Mounting unit
PMDC06



Kit of alluminium
multiadjustable
Hardware with
horizontal bar
LONG PMDC05
SHORT PMDC05C

Hardware



Mounting
examples

Capitis **anatomico**
including
hardware with
long horizontal bar
and detachable
mounting unit


ART.	SIZE
PMDC02XP	Extra Small
PMDC02P	Small
PMDC02M	Medium
PMDC02G	Large

Capitis **anatomico**
including
hardware with
short horizontal
bar and detach-
able mounting unit


ART.	SIZE
PMDC02XPC	Extra Small
PMDC02PC	Small
PMDC02MC	Medium
PMDC02GC	Large

Configurations

Range of depth adjustments




Horizontal bar
upward orientated



Horizontal bar
downward orientated

Accessories



Mounting Interface Plate for backrests

Mounting interface plate for backrests with
plastic or aluminum shell.
The Capitis system can be mounted on solid
backrests only.



Manufactured and distributed by:

 **Promedicare**

Pro Medicare srl
Via A. Montagna - Z.I. - Lotto 41
72023 Mesagne (Br) - Italy
tel. (+39) 0831 77 78 40
info@promedicare.it

www.promedicare.eu

 ProMedicareIt
 _promedicare_
 promedicare



**2 YEARS
WARRANTY**



CERTIFIED COMPANY WITH QUALITY
MANAGEMENT SYSTEM ISO 9001/ ISO 13485

The manufacturer reserves the right to make changes without notice. The contents, the form, the design, the image, the structure of this catalogue are property of Pro Medicare S.r.l. Any improper or unauthorized use is strictly forbidden. Any imitation and/or explicit or not explicit copy can be pursued by the Law.

CERTIFIED
ISO 9001
ISO 13485

