



THE NECESSARY ADDITIONS TO THE POSITIONING

Versa **fixatis**[®]
FASTENERS



Pelvic belt

- ✓ 2 fixing points
- ✓ 2 adjustments
- ✓ Metal front button release



CODE#	SIZE	INDICATIONS
PMD 26S	Small	Paediatric/Teenagers
PMD 26M	Medium	Teenagers/ Adults
PMD 26L	Large	Adults

Provided with individually adjustable pads and padded button protector

Mounting Examples



Technical Information

measurements expressed
in centimeters

CODE#	H	L
PMD 26S	6	17
PMD 26M	7	22
PMD 26L	7	27

H = Height of the pad
W = Width of the pad

Pelvic belt

- ✓ 4 fixing points
- ✓ 2 adjustments
- ✓ Metal front button release



CODE#	SIZE	INDICATIONS
PMD 25S	Small	Paediatric/Teenagers
PMD 25M	Medium	Teenagers/ Adults
PMD 25L	Large	Adults

Provided with individually adjustable pads and padded button protector

Mounting Examples



Technical Information

measurements expressed
in centimeters

CODE#	H	L
PMD 25S	6	17
PMD 25M	7	22
PMD 25L	7	27

H = Height of the pad
W = Width of the pad

Abduction dynamic harness

Promotes the abduction of the legs.

When used in combination with a contoured seat, helps to maintain stable the position of the pelvis. Thanks to the materials it is made, it gives stability allowing at the same time small ranges of movement.




CODE#	SIZE	INDICATIONS
PMD 18S	Small	Paediatric/Teenagers
PMD 18M	Medium	Teenagers/ Adults
PMD 18L	Large	Adults

Mounting examples



Technical Information

measurements expressed
in centimeters

	CODE#	A	B	C
	PMD 18S	28	17	41
	PMD 18M	33	20	49
	PMD 18L	38	24	55

Pelvic belt

- ✓ 2 fixing points
- ✓ 2 adjustments
- ✓ Plastic side squeeze buckle

CODE#	SIZE	INDICATIONS
PMD 19XS	Extra Small	Paediatric
PMD 19S	Small	Teenagers
PMD 19M	Medium	Teenagers/ Adults
PMD 19L	Large	Adults



Provided with individually adjustable pads and padded button protector

Mounting examples



Technical Information

measurements expressed
in centimeters

	CODE#	H	L
	PMD 19XS	4,5	12
	PMD 19S	5,5	17
	PMD 19M	7	22
	PMD 19L	7	27

H = Height of the pad
L = Length of the pad

Midfoot strap padded with hook and loop fixing

Provides support in the mid part of the foot avoiding slipping forward and ensures correct positioning of the foot on the footplate.

The mounting is secured by hook and loop




CODE#	SIZE	INDICATIONS
PMD 15CS	Small	Paediatric
PMD 15CM	Medium	Teenagers/ Adults
PMD 15CL	Large	Adults

Mounting examples



Technical Information

measurements expressed in centimeters



CODE#	A	B
PMD 15CS	21	7
PMD 15CM	27	10
PMD 15CL	30	11

Dynamic butterfly harness

Featuring a soft and dynamic chest stabilizer. It enables the positioning of the trunk and the shoulders and helps small ranges of movements




CODE#	SIZE	INDICATIONS
PMD 16XS	Extra Small	Paediatric
PMD 16S	Small	Teenagers
PMD 16M	Medium	Teenagers/ Adults
PMD 16L	Large	Adults

Mounting examples



Technical Information

measurements expressed in centimeters



CODE#	A	B	C	D	E
PMD 16XS	11	20	17	4	29
PMD 16S	17	23	17	5	34
PMD 16M	20	24	26	5	45
PMD 16L	20	26	30	6	52

WARNING:The performances of the harnesses are carried out when used in combination with a pelvic belt and a suitable back height never higher than that of the shoulders

TRUNK - SHOULDERS
CHEST HARNESS WITH SHOULDERS RETRACTION PADS

Chest harness with shoulders retraction pads

Provides a stable support in the anterior part of chest and optimizes shoulder retraction.
It is multi-adjustable and allows a stable fixing position in the posterior part for a safe use.

CODE#	SIZE	INDICATIONS
PMD 17S	Small	Paediatric
PMD 17M	Medium	Teenagers/ Adults
PMD 17L	Large	Adults




Mounting examples



Technical Information

measurements expressed
in centimeters

	CODE#	L	H
	PMD 17S	16	18
	PMD 17M	18	20
	PMD 17L	23	23

TRUNK - SHOULDERS
TRUNK HARNESS WITH ELASTIC SHOULDERS RETRACTORS
AND CHEST BUCKLE CLOSURE

Trunk harness with elastic shoulders retractors and chest buckle closure

Featuring soft pads and elastic ajustable retractors and allows small ranges of movement.
It can be fastened from the front part or if necessary can be crossed.

CODE#	SIZE	INDICATIONS
PMD 21XS	Extra Small	Paediatric
PMD 21S	Small	Teenagers
PMD 21M	Medium	Teenagers/ Adults
PMD 21L	Large	Adults




Mounting examples



Technical Information

measurements expressed
in centimeters

	CODE#	W	H
	PMD 21XS	4	25
	PMD 21S	4	35
	PMD 21M	6	40
	PMD 21L	6	45



Medical device marked CE
in accordance with the
regulation (UE) 2017/745

Versa
Fixatis®
ELEMENTI DI FISSAGGIO

Manufactured and distributed by:

 **Promedicare**

Pro Medicare srl
Via A. Montagna - Z.I. - Lotto 41
72023 Mesagne (Br) - Italy
tel. (+39) 0831 77 78 40
info@promedicare.it

www.promedicare.eu

 ProMedicareIt
 _promedicare_
 promedicare



**2 YEARS
WARRANTY**



CERTIFIED COMPANY WITH QUALITY
MANAGEMENT SYSTEM ISO 9001/ ISO 13485

The manufacturer reserves the right to make changes without notice. The contents, the form, the design, the image, the structure of this catalogue are property of Pro Medicare S.r.l. Any improper or unauthorized use is strictly forbidden. Any imitation and/or explicit or not explicit copy can be pursued by the Law.

CERTIFIED
ISO 9001
ISO 13485

