

2 A Teschio e denti, 1489.

R.L. 19058 v (B. 41 r; K/P 42 v).

Il vacuo della cassa dell'occhio, e 'l vacuo dell'osso sostenitore della guancia, e quello del naso e della bocca sono d'eguale profondità e terminano sotto il senso comune fiede centrale delle sensazioni per linea perpendicolare. E ciascuna d'esse vacuità ha tanto di profondità quant'è la terza parte del volto dell'omo, cioè dal mento a capegli.



A NEW SOLUTION

Positioning system for the pelvis with modular inserts



Why Pucco ?

- ? FACILITATES TRUNK CONTROL
- ? LIGHTWEIGHT
- ? MULTIFUNCTIONAL



**21 Modular
configurable inserts**



*Free Shape
Padding*



**The alternative to
traditional seating shells**

CLINICAL INFORMATION

- An innovative postural unit for legs-pelvis-trunk reporting system. Completely modular and adjustable to suit many applications, being both easy to handle and extremely lightweight.
- Easy to adjust and adapt to different paediatric therapeutic applications and social lifestyles by encouraging complete integration into daily family life.
- For children with poor postural control and in particular the Hypotonic.
- For very small children who need to be put in seated position.
- For children who have leg length discrepancy, wind-swept, pelvic obliquity, pelvic rotation, posterior and anterior pelvis tilts, and a risk of hip dislocations.

PTS Pelvic Total Support

TECHNICAL FEATURES

- 21 Modular configurable inserts as standard.
- Two layers of foam (flexible open celled foam with low spring back action and reticulate foam for ventilation).
- Rear protective pads (flexible open celled foam with low spring back action).
- Cover.
- Pelvic strap.

SPARE PARTS

- Flat base of construction
- Cover

ACCESSORIES

- Leveling kit for backrest





inserto Pucco


21 configurable modular inserts

- Each postural unit is ordered based upon the size of mobility/wheelchair device

Sizes (cm)	Paediatric			
	UUS short	UUS Tall	US short	US Tall
Width X Length	25x30	25x30	30x38	30x38
Backrest height	16	23	21	28
Seat Width Range	min 13 - max 23	min 13 - max 23	min 20 - max 32	min 20 - max 32
Seat Length Range	min 15 - max 25	min 15 - max 25	min 21 - max 32	min 21 - max 32

- Maximum LOAD (US) : 40 KG

- Available with incontinent and Air Exchange Cover





MANUFACTURED BY:

 **Promedicare**


Via Montagna Z.I. - Lotto 41
72023 Mesagne (BR) - ITALY
Tel. +(39) 0831 777840
Email: info@promedicare.it
www.promedicare.eu

CERTIFIED COMPANY WITH QUALITY
MANAGEMENT SYSTEM ISO 9001 / EN ISO 13485

CERTIFIED
ISO 9001
ISO 13485



DISTRIBUTED BY:



 **2 YEAR WARRANTY**

www.versainserto.com