

Order ☐ Quotation ☐

Invoicing Details	Shipping Details
Company name: _____ Delivery date required: ____ / ____ / ____ Delivery address: _____ P. Code.: _____ Town: _____ Tel.: _____ Fax: _____ # Order: _____ Date of the Order: ____ / ____ / ____	Ship to: _____ _____ CUSTOMER STAMP AND SIGNATURE

Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

Note: Any request for variation in dimensions, materials and type compared to the standard does not have **CE** marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device.


Images are purely illustrative and may not fully reflect reality.

Ermes-Tipco Standard Configuration

- Seat adjustable in width, depth also asymmetrically in adduction and abduction
- Kit of inserts for a customised positioning of the pelvis
- Reclining and height adjustable backrest
- Kit of inserts for a customised positioning of the trunk
- Multi-adjustable and removable containment and positioning kit:
 - pelvis (adjustable in width, depth, height, adduction/abduction)
 - arm support (adjustable in width, depth, height, angle, intra/extra rotation)
- Male interface plate
- Tipco indoor frame adapted with female interface plate:
 - Height adjustable seat by gas spring
 - Tilt in space by gas spring from -5° up to ° +30°
 - 4 swivel castors (2 with locks)
 - Multi-adjustable push handle



Sizes of the user







1 RH	depth of the seat RH	_____
1 LH	depth of the seat LH	_____
2 RH	length of the leg RH	_____
2 LH	length of the leg LH	_____
3 RH	depth of the thorax RH	_____
3 LH	depth of the thorax LH	_____
4	height of the shoulders from the seat plane	_____
5	width of the thorax	_____
6	width of the pelvis	_____

Ermes

- | | | | |
|-------------------------------------|--------------|---------------------------------|------------------------------------|
| ☐ size 2 | color | ☐ blue/grey | ☐ fuchsia/grey |
| ☐ size 3 | color | ☐ blue/grey | ☐ fuchsia/grey |
| ☐ size 4 | color | ☐ wine/grey | |

Note: If *Ermes'* sizes 2/3 are choosen, *Tipco indoor base* will be in white color; if *Ermes'* size 4 is choosen, *Tipco indoor base* will be in grey color

Note: ☐ = optional without up-charge ☐ = optional with up-charge

Type of the seat															
MIDI Seat ☐ lateral/anterior control (standard configuration)								DEEP Seat total control (PTS) <u>optional without</u> <u>up-charge</u>							
☐ Greater closure of the seat-backrest angle ($\pm 5^\circ$)															
Thoracic Supports															
☐ Flat pads						☐ Contoured pads									
		Size (cm)		1	2	3	4			Size (cm)		1	2	3	4
		H		7	7	8	10			H		6	7	8	10
		L		7.5	8.5	10	12			L		9.5	11	13	16
Side	Size	Bracket type*				Bracket shape		Side	Size	Bracket type*				Bracket shape	
☐ RH		☐ Fixed ☐ Swing-Away				☐ A ☐ B		☐ RH		☐ Fixed ☐ Swing-Away				☐ A ☐ B	
☐ LH		☐ Fixed ☐ Swing-Away				☐ A ☐ B		☐ LH		☐ Fixed ☐ Swing-Away				☐ A ☐ B	
*In the technical information at the end of this form you can find the dimensions reachable with the two types of brackets (A type = L/90° Bracket; B type = 2.5 cm Offset Bracket)															
Belts and Harnesses (values expressed in cm, see details in the technical information)															
☐ Chest pad				☐ Dynamic Abduction Harness				☐ Pelvic belt							
 (with thoracic supports only)		Size	HxW				☐ XS (2)* ☐ S (3)* ☐ M (4)*			Size	HxW				
		☐ 1	(8x10)							☐ 1	(6x11)				
		☐ 2	(9x13)							☐ 2	(7x13)				
		☐ 3	(10x13)							☐ 3	(7x15)				
☐ Dynamic Butterfly Harness				☐ Jacket Harness				☐ Contoured Pelvic Belt							
		☐	XS (2)*				☐ XS (2)* ☐ S (3)* ☐ M (4)*			Size	HxW				
		☐	S (3)*							☐ 1	(6x11)				
		☐	M (4)*							☐ 2	(8x15)				
										Value H is referred to the height of the central part of the pad					
*The number in the brackets is referred to the relevant Ermes size. Therefore it is advisable to combine the default size.															
Hardrest with multi-axial joint and multi-adjustable arm															
Confort with ☐ detachable hardware ☐ fixed hardware				Anatomic with ☐ detachable hardware ☐ fixed hardware				3 piece pads supportive headrest with ☐ detachable hardware ☐ fixed hardware							
		☐	S				☐ XS ☐ S ☐ M			☐ 1	(7x13)				
		☐	M							☐ 2	(9x20)				
		☐	G							values expressed in cm of the occipital pad only					
Other accessories															
		Knee elevating legrest ☐ RH ☐ LH (with cover for angle knee)													
		separate footplates: ☐ RH ☐ LH													
		Footplate options: slotted separate footplates: ☐ RH ☐ LH													
		☐ footboard (it does not allow asymmetry, adduction and abduction)													


☐ Shoes Holders




☐ Calf support ☐ RH ☐ LH
☐ Calf panel (it does not allow asymmetry, adduction and abduction)
☐ Lateral supporting pad for legs ☐ RH ☐ LH
☐ Padded anterior/posterior Calf strap
☐ Padding for separate footplates
☐ Padding for footboard


☐ Elbow lateral-posterior support ☐ RH ☐ LH

☐ Grey detachable Tray
with body contouring recess and edges
☐ standard version
☐ larger version for daily activities

☐ Transparent detachable Tray
with body contouring recess and edges
(only standard version)

Note: the dimensions of the tray in the standard grey color version and those of the transparent tray are the same, the naming changes (the size 2 of the grey tray will be the size 1 of the transparent tray, 3->2, 4->3)

standard version dimensions				
	Size (cm)	L	P	P1
2 (for Ermes size 2)	40	35	12	
3 (for Ermes size 3)	48	40	13.5	
4 (for Ermes size 4)	58	45	15.5	

larger version for daily activities dimensions				
	Size (cm)	L	P	P1
1 (for Ermes size 2)	50	45	13,7	
2 (for Ermes size 3)	60	50	15.4	
3 (for Ermes size 4)	68	55.6	19.3	

☐ Padding of the tray ☐ upper ☐ lower

☐ **Additional Outdoor Frame Adacta Giò for Ermes**

Standard Configuration

- Adacta titanium satin anodised frame which plates have the same color of the cover
- Adjustable single handlebar with folding backrest tubes
- Reclining mechanism by predetermined plate w/pin 3°, 10°, 18°
- Continuous Tilt in space by drive lever: from 2° up to 35° (from 6° up to 39°) tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks
- Fixed anti-tip system
- Wheels Kit with brake
- Height adjustment: high (49 cm from the floor)/low (43 cm from the floor) with 300mm and 500mm wheels only
- Solid grey base assembled with female interface plate





◇ = Standard **○ = Optional without up-charge** **□ = Optional with up-charge**

Frame size	2 (combined with Ermes size 2)	3 (combined with Ermes size 3)	6 (combined with Ermes size 4)
Castors	◇ 150mm (6") poly grey ☐ black	☐ 175mm (7") poly	☐ grey ☐ black
Ruote posteriori	◇ 300mm (12") poly grey with parking brake lever ☐ 300mm (12") poly black with parking brake lever ☐ 300mm (12") pneumatic grey with parking brake lever		
	☐ 300mm (12") poly black with drum brake ☐ 300mm (12") poly grey with drum brake ☐ 300mm (12") pneumatic grey with drum brake		
	☐ 300mm (12") poly black with drum brake and parking brake lever ☐ 300mm (12") poly grey with drum brake and parking brake lever ☐ 300mm (12") pneumatic grey with drum brake and parking brake lever		
	☐ 400mm (16") poly black with parking brake lever (only 45 cm height) ☐ 400mm (16") poly black with drum brake (only 45 cm height) ☐ 400mm (16") poly black with drum brake and parking brake lever (only 45 cm height)		
	☐ 500mm (20") poly grey without handdrims with parking brake lever ☐ 500mm (20") pneumatic grey without handdrims with parking brake lever ☐ handdrims for 500mm (20") wheels		
	☐ 500mm (20") pneumatic grey with handdrims with drum brake ☐ 500mm (20") pneumatic grey with handdrims with drum brake and parking brake lever		

Back Posts	☐ Kit for fixed push handles ☐ Kit for telescopic and angle adjustable single handlebar ☐ Kit for telescopic push handles with lever block (max adjustment 10 cm) ☐ Kit for angle adjustable separated handlebar ☐ Quick-release pins for easier back reclining
Tilt in space	☒ <i>with pedal</i> ☐ <i>with lever</i>
☐ Canopy grey/black	adaptor included
☐ Basket	
☐ Rain cover	
Foot muff	☐ Blue ☐ Pink ☐ Bordeaux

TECHNICAL INFORMATION

Ermes		Size 2	Size 3	Size 4
Width of the Seat (cm)		≤25	18-28	22-35
Depth of the Seat (cm)		19-26	22-34	31-43
Height of the trunk (cm)		26-38	33-45	41-53
Width of the thorax (cm)		bracket A 19-25	bracket A 22-28	bracket A 27-33
		bracket B 14-20	bracket B 17-23	bracket B 22-28
Max abduction of the leg (°)		16	16	16
Max adduction of the leg (°)		16	16	16
Popliteal Space to footplate Length (cm)* <i>*if the elevating knee angle is provided</i>		8-19	8-29	8-39
Height of the arm support (cm)		14-19	20-26	20-26
Backrest Reclining (°)		-3+20	-3+20	-3+20
Ermes load (kg)		25	35	50
Ermes weight (kg)		5.2	8.1	12.1
Contoured hip guide		Size 2	Size 3	Size 4
	H	6.5	7	8
	P	17	20	23
Flat arm support		Size 2	Size 3	Size 4
	L	7.5	7.5	8.5
	P	17	20	23
Note: the postural unit for armrest and the lateral hip pad are combined to the relative size of Ermes and are not interchangeable				
Tipco		Size 2	Size 2	Size 3
Seat to floor height (cm)		63	63	63
Width overall dimension (cm)		67	67	67
Depth overall dimensions (cm)		80	80	88
Tilt in space (°)		-5+30	-5+30	-5+30
Tipco Load (kg)		50	50	50
Tipco Weight (kg)		13.6	13.6	14
Total Weight (Tipco+Ermes) (kg)		18.8	21.7	26.1
Adacta Giò		Size 2	Size 3	Size 6
Seat to floor height (cm)		43	43	43
Width overall dimension (400/500 wheels) (cm)		48	53	58
Min/max depth overall dimensions (cm)		73-100	75-102	81-109
Tilt in space (°)		2-35 (from 6 up to 39) <i>tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks</i>		
Predetermined backposts reclining (°)		3-10-18		
Adacta Giò load (standard configuration) (kg)		30	50	50
Adacta Giò weight (kg)		12.8	13.4	14.5
Total weight (Adacta Giò+Ermes) (kg)		18	21.5	26.6