



Order ☐ Quotation ☐






Invoicing Details	Shipping Details
Company name: _____ Delivery date required: ____ / ____ / ____ Delivery address: _____ P. Code.: _____ Town: _____ Tel.: _____ Fax: _____ # Order: _____ Date of the Order: ____ / ____ / ____	Ship to: _____ _____ CUSTOMER STAMP AND SIGNATURE



Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

Note: Any request for variation in dimensions, materials and type compared to the standard does not have **CE** marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device.

Images are purely illustrative and may not fully reflect reality.

Standard Configuration																									
- Adacta titanium satin anodised frame with colored plates - Adjustable single handlebar with folding backrest tubes - Reclining mechanism by predetermined plate w/pin 3°, 10°, 18° - Continuous Tilt in space by pedal: from 2° up to 35° (from 6° up to 39°)¹ ¹ tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks - Fixed anti-tip system - Wheels Kit with brake - Adjustment of the seat depth - Height adjustment: high (49 cm from the floor)/ low (43 cm from the floor)* *(with 300mm and 500mm wheels only) - Knee Angle - Foot Plates - Simple Calf Strap																									
◆ = Standard ○ = Optional without up-charge □ = Optional with up-charge Note: If a selection is not made, the product will be delivered in its standard configuration. (See bold items)																									
Rigid Frame Structure																									
Frame Color	Titanium satin anodised frame with colored plates: ○ Blue ○ Pink ○ Bordeaux																								
Frame size (Width measurements are from tube to tube externally, while depth measurements are from back post to footboard engagement)																									
<table border="1"> <thead> <tr> <th>Size</th> <th>Seat Width (cm)</th> <th>Depth and adjustments (cm)</th> </tr> </thead> <tbody> <tr> <td>○ 1</td> <td>25</td> <td>○ 32 ◆ 34.5 ○ 37</td> </tr> <tr> <td>○ 2</td> <td>25</td> <td>○ 36 ◆ 38.5 ○ 41</td> </tr> <tr> <td>○ 3</td> <td>30</td> <td>○ 32 ◆ 34.5 ○ 37</td> </tr> <tr> <td>○ 4</td> <td>30</td> <td>○ 36 ◆ 38.5 ○ 41</td> </tr> <tr> <td>○ 5</td> <td>35</td> <td>○ 36 ◆ 38.5 ○ 41</td> </tr> <tr> <td>○ 6</td> <td>35</td> <td>○ 40 ◆ 42.5 ○ 45</td> </tr> <tr> <td>○ 7</td> <td>38</td> <td>○ 40 ◆ 42.5 ○ 45</td> </tr> </tbody> </table>	Size	Seat Width (cm)	Depth and adjustments (cm)	○ 1	25	○ 32 ◆ 34.5 ○ 37	○ 2	25	○ 36 ◆ 38.5 ○ 41	○ 3	30	○ 32 ◆ 34.5 ○ 37	○ 4	30	○ 36 ◆ 38.5 ○ 41	○ 5	35	○ 36 ◆ 38.5 ○ 41	○ 6	35	○ 40 ◆ 42.5 ○ 45	○ 7	38	○ 40 ◆ 42.5 ○ 45	
Size	Seat Width (cm)	Depth and adjustments (cm)																							
○ 1	25	○ 32 ◆ 34.5 ○ 37																							
○ 2	25	○ 36 ◆ 38.5 ○ 41																							
○ 3	30	○ 32 ◆ 34.5 ○ 37																							
○ 4	30	○ 36 ◆ 38.5 ○ 41																							
○ 5	35	○ 36 ◆ 38.5 ○ 41																							
○ 6	35	○ 40 ◆ 42.5 ○ 45																							
○ 7	38	○ 40 ◆ 42.5 ○ 45																							
	H: seat rails to floor height (see page 4) Q: length from seat rails to footplate (see page 4) C: push handles height to floor = 110 cm M: push handles height from seat rails (see page 2)																								

Vehicle Transportation Kit (according to the ISO 7176-19 standard) Do not use the wheelchair as a seat in a vehicle if the weight of the user is lower than 22 kg ☐ only frame height 43 cm and with 300mm (12") and 500mm (20") posterior wheels <i>(it is strongly recommended the use of a headrest)</i>									
Wooden Supporting Base for Breathing Equipment <i>available with tilt in space with lever only</i> Note: Bulky equipments may lead to a reduction of the tilt in space /reclining values ☐ covered in black leather flex Please report sizes: Width _____; Length _____; Height _____; NOTE _____ Note: Pro Medicare Srl reserves the right to check wheelchair compatibility									
Pedal for barriers									
Side ☐ RH ☐ LH									
Leg Support									
Separate Leg Cradle ☐ RH ☐ LH <i>adjustable in height and depth</i>									
	cm	1	2	3	4	5	6	7	
	H	6	6	10	10	12,5	12,5	15	
	L	6	6	9,5	9,5	12,5	12,5	12,5	
☐ Calf Panel <i>adjustable in height and depth</i>									
	cm	1	2	3	4	5	6	7	
	H	8	8	10	10	12	12	15	
	L	17	17	23	23	26,5	26,5	25	
☐ Padded calf strap									
	cm	1	2	3	4	5	6	7	
	H	11	11	11	11	11	11	11	
	L	20	20	26	26	31	31	34	
Cover for Knee Angle joint ☐ RH ☐ LH					Cover for legrest tube ☐ RH ☐ LH				
Armrest									
☐ Height adjustable flip back tubular armrests/inclination max 5° (L = 25 cm)									
☐ Kit of ABS laterals multi-adjustable, detachable and w/multi-adjustable arm support <i>in width and depth with 2.5 cm gap</i> <u>make a selection:</u> ☐ small arm support (from 1 to 5) (dimensions: 5 cm x 25 cm) ☐ large arm support (from 6 to 7) (dimensions: 6 cm x 30 cm) <u>make a selection:</u> ☐ short vertical bar (from 22 cm to 31 cm*) ☐ long vertical bar (from 29 cm to 38 cm*) * distance from armrest to the seat rail									
☐ Multi-adjustable and detachable hips and arms positioning kit - pelvis: width, depth, height, adduction/abduction adjustable - arm support: width, depth, height, inclination, intra/extra-rotation adjustable <i>For independent adjustments referring to the width the gap is 3 cm</i> The armrests and the hip guides are assembled on the frames as the following indications without possibilities to be modified ☐ small (from frame size 1 up to 4) (armrests dimensions: 8 cm x 22 cm) (hips guide dimensions: 25 cm x 10 cm with slot of 5 cm) ☐ large (from frame size 5 up to 7) (armrests dimensions: 10 cm x 30 cm) (hips guide dimensions: 30 cm x 12 cm with slot of 7.5 cm) 									
									
Tray (recommendation for choice: size 1 for Giò size 1 and 2, size 2 for Giò size 3 and 4, size 3 for Giò size 5, 6 and 7)									
grey wooden with body contouring recess and edges ☐ standard option ☐ larger option for daily activities					with Hardware ☐ w/lateral flip ☐ detachable				
☐ transparent with body contouring recess and edges <i>(only standard option)</i>					with Hardware ☐ w/lateral flip ☐ detachable				
Padding ☐ upper ☐ lower					 lateral flip		 detachable		

standard option dimensions						larger option for daily activities dimensions					
	Size (cm)	P	L	B	A		Size (cm)	P	L	B	A
	1	40	48	28	13		1	45	52.5	28	13.5
	2	45	58	35	15		2	50	62.5	36	15.5
	3	50	66	40	19		3	55	71	42	19.5

VERSA-FIXATIS FASTENERS

☐ **Pelvic Belt** with 2 fixing points and 2 adjustments with plastic side squeeze buckle

	cm	☐ XS	☐ S
	H (pad height)	4.5	5.5
	L (pad length)	12	17

☐ **Pelvic Belt** with 2 fixing points and 2 adjustments with metal push button buckle

	cm	☐ S	☐ M
	H (pad height)	6	7
	L (pad length)	17	22

☐ **Pelvic Belt** with 4 fixing points and 2 adjustments with metal push button buckle

	cm	☐ S	☐ M
	H (pad height)	6	7
	L (pad length)	17	22

OTHER ACCESSORIES

 ☐ Canopy grey/black adaptor included	 Foot Muff ☐ Blue ☐ Pink ☐ Bordeaux	 ☐ Rain cover	 ☐ Basket
--	--	---	--

FRAME DIMENSION AND FEATURES

Size	1	2	3	4	5	6	7
Load (kg)/Frame Weight* (kg)	30/12.6	30/12.8	50/13.4	50/13.6	50/14.1	50/14.5	50/15.1

* Weight of the frame, including wheels, knee angle and footplates as standard configuration

Weight and overall dimensions (cm) of the frame with backrest tubes folded down, footrest and posterior wheels detached

Size	1	2	3	4	5	6	7
Length	64	66	64	66	66	68	68
Width	38	38	42	42	48	48	51
Height	64	64	65	65	66	66	67
Weight (kg)	9.4	9.6	10.1	10.3	10.7	11.1	11.5

Note: The length dimensions, referred to the average values, are subjected to the variations due to the displacement of the wheel plates. The values reported are intended for the standard configuration

Width (cm)

Size	1	2	3	4	5	6	7
300mm wheel	standard	43	43	48	48	53	56
	w/drum brake	44	44	49	49	54	57
400/500mm wheel	standard	45	45	50	50	55	58
	w/drum brake	48	48	53	53	58	61

