

leo

VIP



Comfort and **prevention** for the user,
convenience of use for the caregiver.

promedicare.eu

Risposte reali a bisogni speciali
ProMedicare®



The all-in-one solution designed to maintain the correct and safe posture. It combines innovative technologies with an ergonomic design, ensuring containment, comfort and prevention for the user, convenience of use and lightness for the caregiver.

The all-in-one comfort solution



Headrest with multi-adjustable bracket



Tilt and height adjustable single handlebar



Pelvic belt



With Tensioning strap system



Waterproof, antibacterial, fireproof and breathable double cover



- Loads distribution
- Lateral containment
- Ischial, trochanteric, coccygeal and sacral suspension

- Without mounting hardware and aluminium shell
- Maximum support, control and containment of the trunk

Swing away and detachable knee angle



Standard Configuration

Outdoor Frame for Positioning System

- Leo outdoor frame for positioning system reclining and tilting with gas piston and 43 cm height depth adjustable
- Height, tilt-adjustable, folded down backrest tubes with tilt and height adjustable single handlebar
- Anti-tip wheels
- Parking brake lever
- 175 mm (6") poly casters black
- 300 mm (12") posterior wheels black
- Detachable arm support adjustable in height, width and depth for ABS laterals
- 75° swing-away and detachable knee angle
- Fixed wooden seat base with slots for seat depth adjustment

Positioning Systems

- Inserto Seat Vip pelvis positioning system
- Inserto Back Vip trunk positioning system with tensioning strap system
- Capitis Confort head positioning system size M with multi-adjustable bracket
- Pelvic belt with plastic buckle size M

For the complete list of the optional parts and/or accessories, please refer to the latest order form in force.

		Frame Dimensions				
		XS1	S	M	M1	M1*
Frame	Length (cm)**	71	73	83	83	83
	Width (cm)**	51	54	54	57	57
	Height (cm)**	67	67	68	68	68
	Weight (Kg)**	14,5	15,1	15,7	16,7	16,7
	Load	75	75	100	100	100
	Frame Weight (Kg)***	19,8	20,4	21	21,7	21,7
Cushion (width x depth) cm		40x40 40x45	42x42 42x45	42x45 42x50	45x45 45x50	48x48 48x50
Backrest (width x height) cm		40x45	42x45	42x45 42x50	45x45 45x50	48x50

* with 1.5 cm spacers on each side

** Weight and dimensions of the frame with backrest tubes (fixed) folded down, footrest and posterior wheels detached

*** Weight of the frame, including wheels, knee angle and footplates as per standard configuration



Manufactured and distributed by:



Pro Medicare srl
Via A. Montagna - Z.I. - Lotto 41
72023 Mesagne (Br) - Italy
tel. (+39) 0831 77 78 40
fax (+39) 0831 73 07 39

info@promedicare.it
www.promedicare.eu
f ProMedicareIt
i promedicare_
v promedicare



**2 YEARS
WARRANTY**



CERTIFIED COMPANY WITH QUALITY
MANAGEMENT SYSTEM ISO 9001/ ISO 13485

The manufacturer reserves the right to make changes without notice.
The contents, the form, the design, the image, the structure of this catalogue are property of Pro Medicare S.r.l. Any improper or unauthorized use is strictly forbidden. Any imitation and/or explicit or not explicit copy can be pursued by the Law.

